

TECHNICAL FEATURES

- Prismatic telescopes with magnification from 3x to 10x
- patented interpupillary distance adjustment system
- angle of inclination of the telescopes adjustable from 12 ° to 35 ° with respect to the visual axis
- possibility of inserting corrective lenses in telescopes and frames for additional lenses
- optical unit made and collimated at the required working distance, between 30 and 60 cm
- headband with six rotative points for better adaptation to the skull shape.
- optical group regulated by three independent micrometric movements that allow the overturning on the front without losing the chosen position.

COMPATIBLE ACCESSORIES

- LED light with rechargeable battery
- removable corrective lenses
- removal frame for corrective lenses

to obtain tri-focal system: far view,
intermediate view, close view

- magnetic system to fix LED light
- protective flip up mask