


SLH Prismatic Loupes




- ▶ Pro Headbands can be mounted a cordless battery pack. User can choose either rear band or top band mounted.
- ▶ Achromatic lenses give a clear image edge to edge with no aberration. Choice of 4x, 5x or 6x magnification.
- ▶ Independently adjustable optics ensure accurate setting. Choice of working distance 340,420,460,500mm.
- ▶ Coated optics minimise reflections and provide crisp, bright images.
- ▶ Optics are fully adjustable for height and declination.
- ▶ Water resistant for easy cleaning.
- ▶ Optional LED illumination with cordless battery pack.



Locker with one key



Interface with magnet



Sterilizing cover

Magnification	Working distance	Field of view	Magnification	Working distance	Field of view
4.0x	340mm/13"	55mm/2.2"	5.0x	340mm/13"	46mm/1.8"
4.0x	420mm/17"	65mm/2.6"	5.0x	420mm/17"	55mm/2.2"
4.0x	460mm/18"	75mm/3.0"	5.0x	460mm/18"	60mm/2.4"
4.0x	500mm/20"	85mm/3.3"	5.0x	500mm/20"	65mm/2.6"

Magnification	Working distance	Field of view
6.0x	340mm/13"	40mm/1.6"
6.0x	420mm/17"	47mm/1.8"
6.0x	460mm/18"	50mm/2.0"
6.0x	500mm/20"	58mm/2.3"

Optional 3W HL8200 LED light

- 3W LED loupe light, homogenous well defined illumination. Approx. 55,000 Lux.
- Quick recharge time of less than 2 1/2 hours. Up 6 hours continuous use.
- Removable and washable comfort pads, front, rear and top.
- With a cordless battery pack.



SLH Prismatic Loupes with professional headband, Cordless HL8200 headlight, Lithium ion Headband mounted battery pack, charging dock, power cable, transformer and carry case. Battery pack can be top mounted or rear mounted.



Cordless HL8200 headlight, Lithium ion Headband mounted battery pack, charging dock, power cable, transformer and carry case. Battery pack can be top mounted or rear mounted.



Spare lithium ion headband mounted battery pack.



Orange filter for use in Dentistry to help prevent premature composite hardening.



Sterilizing cover