

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +39 0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



TECHNICAL SHEET

PRODUCT NAME	PROTECTIVE EYESHIELDS WITH AUTOCLAVABLE FRAME
CODE	R66
COMPONENT	PROTECTIVE VISOR: POLYCARBONATE (SABIC) VISOR THICKNESS: 0.250 mm FRAME: MAKROLON 2458 (POLYCARBONATE)
THERMAL CHARACTERISTIC	FRAME AUTOCLAVABLE UP TO 121°C.
PHYSICAL CHARACTERISTICS, MECHANICAL CHARACTERISTICS AND PACKING	PACKING IN CARTON BOX: 1 FRAME AND 6 PROTECTIVE SHIELDS. THE VISOR COVER NOSE MOUTH AND EYES. THE PROTECTIVE VISOR HAS A DOUBLE PROTECTIVE COATING. THE VISOR IS VERY TRANSPARENT.
CHEMICAL CHARACTERISTIC	IT RESISTS TO THE CHEMICAL DISINFECTANTS.
SAFETY DATA SHEET	N/A
PRECAUTION	WARNING: USE ONLY AS FACE SHIELD AGAINST SPLASHING LIQUID OF WATER SALIVA.
INSTRUCTIONS FOR USE	SUPPLIED IN "NON-STERILE" CONDITION. AUTOCLAVABLE. AVAILABLE INSTRUCTION FOR FITTING INSIDE CARTON BOX.
ANALYSIS OF THE RISKS REF.D.L.46 24/2/97	AVAILABLE ONLY IN ITALIAN
CLAS. CEE EN 166	IS AN INDIVIDUAL PROTECTION DEVICE(DPI) CLASS I

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +39 0185 938 595

info@larident.it - www.larident.it

LARIDENT

MADE IN ITALY
SINCE 1985



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	VISIERA PROTETTIVA IN POLICARBONATO CON MONTATURA AD INCASTRO
NS CODICE	R66
MATERIALE	VISIERA: POLICARBONATO (SABIC) SPESSORE VISIERA: 0.250 mm MONTATURA: MAKROLON 2458 (POLICARBONATO)
PROPRIETÀ TERMICHE	LA MONTATURA E' STERILIZZABILE A 121°C.
PROPRIETÀ MECCANICHE, FISICHE, CARATTERISTICHE, CONFEZIONAMENTO	LA MONTATURA E' DISPONIBILE IN 5 COLORI. LA CONFEZIONE: 1 MONTATURA+6 VISIERE CHE COPRONO NASO OCCHI E BOCCA, HA UN DOPPIO RIVESTIMENTO PROTETTIVO ED E' DOTATA DI UN ECCELLENTE TRASPARENZA.
PROPRIETÀ CHIMICHE	RESISTE A TUTTI GLI AGENTI CHIMICI DISINFETTANTI.
SCHEDA DI SICUREZZA	N/A
PRECAUZIONI	ATTENZIONE: GLI OCCHIALI PROTETTIVI PROTEGGONO SOLO ED ESCLUSIVAMENTE DA GOCCE E SPRUZZI D'ACQUA O SALIVA.
ISTRUZIONI PER L'UTILIZZO	DESTINATI AD ESSERE UTILIZZATI ESCLUSIVAMENTE DAL MEDICO SPECIALISTA O DA PERSONALE ALTAMENTE SPECIALIZZATO. VIENE FORNITO IN CONDIZIONE "NON STERILE".
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97, 37/2010	SI
CLASSIFICAZIONE CEE EN 166	DISPOSITIVO DI PROTEZIONE INDIVIDUALE (DPI) CLASSE I

N/A = NON APPLICABILE

N/N = NON NECESSARIA